

BOECO AUTOCLAVE BTE-23D

Features:

- ▶ Three times fractionated pre-vacuum
- ▶ Stainless steel chamber
- ▶ LCD display
- ▶ Equipped with high quality and high performance pump
- ▶ 1 tray holder, 1 tray handle and 3 trays included
- ▶ High quality water pump and valves
- ▶ Independent steam generator
- ▶ Two samples self-calibration
- ▶ Overhead type water storage tank, easy to clean.
- ▶ Drying procedure: Dry by vacuum.
- ▶ Sterilization temperature: 121°C and 134°C
- ▶ Test programs: Helix test, B&D test and Vacuum test.
- ▶ Sterilization cycles: 7 cycles are available and 3 test programs.

Program	Temperature	Pressure	Sterilization time	Total time
SOLID 134	134°C	2.1 bar	4 min	18-30 min
SOLID 121	121°C	1.1 bar	20 min	30-45 min
TEXTILE 134	134°C	2.1 bar	6 min	30-40 min
(or wrapped)				
TEXTILE 121	121 °C	1.1 bar	20 min	35-50 min
(or wrapped)				
LIQUID 134	134°C	2.1 bar	10 min	45-65 min
LIQUID 121	121°C	1.1 bar	30 min	50-75 min
PRION	134°C	2.1 bar	18 min	45-65 min
B&D TEST	134°C	2.1 bar	3.5 min	22-35 min
Helix TEST	134°C	2.1 bar	3.5 min	22-35 min
Vacuum TEST	-	-	-	15-20 min

- ▶ Only have several buttons, very easy to use, any people can run it easily.
- ▶ All cycles are fully automatic; there is no need for human intervention in cycle management.
- ▶ Pressure protection locking system that ensures the door will not be opened once there is any pressure in chamber.
- ▶ Sterilization documentation: Integrated Micro Printer and USB port to store the information of the programm cycle in flash memory.
- ▶ The sterilization time and dry time can be adjusted by manual, ensure user obtain best sterilization effects according to private requirements.
- ▶ Malfunction detection system, corresponding error code will be showed on the display once there's something goes wrong.

Safety:

- ▶ Warning system prevents starting a cycle if door is not properly locked.
- ▶ Triple Threat Protection for overheating and overpressure.
 - The pressure protection system will automatically cut off the electricity supply once the temperature of chamber inside is abnormally. Motorola pressure sensor is used to accurately sense the state of pressure inside chamber.
 - Temperature control system is using high quality PT1000 temperature sensor, it is high precision, reliable service and prompt measure, ensure accurate and timely transmission of temperature fluctuation.
 - The safety valve will release steam in case the pressure inside chamber exceed normal levels.


Technical Information:

- ▶ Conform to CE and EN13060
- ▶ Manufactured under ISO 13485 norms
- ▶ Chamber Size: ø 247 x 460 mm (23 L)
- ▶ External dimension: 690 x 480 x 455 mm, Netweight 51 kg
- ▶ Packing box 750x560x460 mm, Grossweight 58 kg
- ▶ Power: 220V, 50/60 Hz or 110V, 60HZ

Code	Description
BOE 8902300	BOECO BTE-23D Bench Top Autoclave, Vol. 23 L. incl 1 tray, 1 tray handle and 3 trays, 220V, 50/60Hz
BOE 8902310	BOECO BTE-23D Bench Top Autoclave, Vol. 23 L. incl 1 tray, 1 tray handle and 3 trays, 110 V, 60 Hz



BOECO BTE-23D



MICRO PRINTER AND USB PORT

This BTE-23D adopts 3 times pre-fractionated vacuum system, the following chart simply describe the status of its process

