

# FMH laboratory products

QUALITY  
GUARANTEED

## Cooled Baths

- \*Temperatures down to - 20°C \*
- \*Safety thermostat \*
- \*Drain valve \*
- \*Two year warranty \*

FMH Cooled water baths conform to the Low Voltage Directive 73/23/EEC and are manufactured to the highest standards with the following features:-

- \* High grade stainless steel tank and top \*
- \* Stainless steel flat lid \*
- \* Drain valve at rear of unit\*
- \* Over temperature safety thermostat with indicator\*
- \* Front panel on/off switch\*
- \* Robust epoxy coated housing\*
- \* Recessed front control panel\*
- \* Two year warranty\*
- \* Comprehensive operating manual\*

Options:-

- 1) -40°C units

Note : Circulating temperature controller optional. (see TRE brochure for specifications )

Special requirements on request.



TNK14R20  
+ TR-E / L

MODEL	TNK14R20	TNK23R20
Temperature - Min	-20°C	-20°C
Tank volume (Nominal) liter	14	23
Tank volume (Brim) liter	18	30
Cooling Capacity @ 0°C (Watt)	300	300
Tank aperture	300x300	500x300
Aperture with FMH Controller	300x165	300x365
Dimensions (mm) LxWxH	410x480x530	660x480x530
Nett weight (kg)		
Power supply 50Hz	230V±10%	230V±10%

Specifications :- ambient @ 25°C.

Due to the continuous development of our products, the goods supplied may vary in detail to that illustrated in this brochure.



FMH electronics cc  
P.O.BOX 5187  
HELDERBERG 7135  
CELL. 0829225934  
E-MAIL, fmh7@global.co.za  
TEL - 021 8555330  
FAX - 021 8555484

CK1999/037262/23

Supplied by:-