

# FMH laboratory products



## Tanks

- \* Insulated \*
- \* Two year warranty\*

FMH Insulated Tanks are manufactured to the highest standards with the following features:-

- \* High grade 304 stainless steel inner tank
- \* Epoxy coated housing \*
- \* Insulated \*



TNK24 with Lid

### Options

- \* Height reducers
- \* Adjustable Height reducers \*
- \* Gabled and Flat Lids \*



Tank with controller  
(circulating)

MODEL	Depth Max (mm)	Nominal Capacity (l)	Brim Capacity (l)	Aperture (mm) LxW	Aperture with bridge plate LxW	Dimensions External (mm) LxWxH	Weight Nett Kg
TNK 15	200	13	18	300x300	165x300	385x360x240	
TNK 24	200	18	24	400x300	265x300	485x360x240	
TNK 34	200	25	35	550x300	415x300	635x360x240	
TNK 40	200	30	40	710x300	575x300	795x360x240	

Specifications :- ambient @ 25°C.

**Special Tanks, Refrigerated tanks, Lids etc- available.**

Due to the continuous development of our products, the goods supplied may vary in detail to that illustrated in this brochure.



**FMH** electronics cc  
 P.O.BOX 5187  
 HELDERBERG 7135  
 CELL. 0829225934  
 E-MAIL, fmh7@global.co.za  
 TEL - 021 8555330  
 FAX - 021 8555484

CK1999/037262/23

Supplied by:-