



The Scientific range of laboratory steambath or hot water baths are ideally suited for many types of digestion experiments. The self locating removable lids can readily accommodate a variety of round bottom kjeldahl flasks. All models are manufactured from high grade 304 stainless steel. The unique design of the lid and tank rim ensure that all water vapour returns back to the bath.



Model 402



Model 401

Features:

- Easily replaceable high quality incolloy heating element for long and reliable use
- Manual reset safety thermostat protects the element should the bath run dry
- Simple device for automatic control of water level
- Flask openings 31mm, 50mm, 70mm, and 80mm diameter
- Heat resistant, safe to touch lid knobs

Technical data - Steambaths

Model No.	401	402
Number of stations.	6	12
External dimensions (mm)		
Width	410	530
Height	110	110
Depth	280	405
Electrical connection (V/hz)	230/50	230/50
Power rating @ 230V (W)	2000	2000
Gross weight (kg)	5	7