

# Ovens

- \* Micro Processor Models \*
- \* Thermostatic controlled Models\*
- \* Two year warranty\*

Ovens have 304 stainless steel inners, powder coated outer housings and fitted with a temperature controller ,digital or thermostatic plus the following:- over temperature safety thermostat \* over temperature indicator lamp \* fan assisted air circulation \* non tip self runners \* adjustable fresh air vent and supplied with two chrome plated wire shelves

**Options:-** \* Perforated S/S shelves



ODS 80

## OVENS, THERMOSTAT - ECONOMY

Model	Volume (Litre)
OTE 40	40
OTE 80	80
OTE 160	160
OTE 240	240
OTE 400	400

Temperature Control range : ambient +5°C to 250°C (min 30°C)  
 Control accuracy : +/- 2°C at 150°C  
 Chamber uniformity : 1°C at 100°C  
 Specifications :- ambient @ 25°C.

**Options:-** Digital timer - Max 9999H Switches heater off

## OVENS, DIGITAL - STANDARD

Model	Volume (Litre)
ODS 40	40
ODS 80	80
ODS 160	160
ODS 240	240
ODS 400	400

Temperature Control range : ambient +5°C to 250°C  
 Control accuracy : +/- 0,2°C at 100°C  
 Chamber uniformity : 0,5°C at 100°C  
 Specifications :- ambient @ 25°C.

**Options:-** Digital timer - Max 9999H Switches heater off  
 Ramp controller patterns x segments - max 64

MODEL Capacity	No. of Shelf Max	Dimensions Internal (mm) HxWxD	Dimensions External (mm) HxWxD	Wattage (Kw)	Weight (Kg) Approx.
40 L	3	350x350x330	470x570x540	1,0	43
80 L	4	450x450x400	570x670x610	1,5	60
160 L	5	550x580x500	670x800x710	2,0	95
240 L	6	650x720x520	770x940x730	2,5	118
400 L	10	1050x720x520	1170x940x730	3,0	153

Due to the continuous development of our products, the goods supplied may vary in detail to that illustrated in this brochure.