

Heating Block Units

- * Temperature range... 0 to 450°C *
- * Micro Processor PID Controller *
- * Digital Temperature indicator *
- * Audio / Visual High Alarm *
- * Timer 1 to 9999 minutes *
- * Ramp controller model *
- * Two year warranty *

FMH Dry Block Heaters conform to the Low Voltage Directive 73/23/EEC and are manufactured to the highest standards with the following features:-

Micro Processor Models Standard

- * Micro Processor PID Controller *
- * Digital temperature indicator *
- * Audio / Visual high temperature alarm *
- * Temperature Deviation indicators
- * Digital Timer and buzzer 1 to 9999 minutes


Micro Processor Models with Ramp

- * Eight segment two pattern Controller *

All Models

- * High grade 6082 machined aluminium blocks - Hole depth 60mm*
- * Stainless steel sheathed heating element*
- * Slanted front control panel*
- * Front panel on/off switch*

Note:- Select a Control unit and a Block unit for a complete system.



MODEL Controller	ABH -TC (standard)	ABH -TC/R (ramp)
Temperature control range(Ambient +5)	0 to 450°C	0 to 450°C
Control accuracy @ 100°C	+/- 0,2°C	+/- 0,2°C

MODEL Block unit	Hole dia (mm)	No. of Holes
ABH-BU	16,5	105
ABH-BU	23	60
ABH-BU	28	40
ABH-BU	47,5	18

Specials on request- beakers, slide etc
Specifications :- ambient @ 25°C.

MODEL	Description	Dimensions External (mm)	Wattage Kw	Weight Nett Kg
ABH-TC	Standard MP Control Unit	175 x 190 x 110H	0,1	
ABH-TC/R	Ramp Control Unit	175 x 190 x 110H	0,1	
ABH-BU	Block unit (Specify hole Dia)	565 x 270 x 165H	2,0	

Due to the continuous development of our products, the goods supplied may vary in detail to that illustrated in this brochure.