

## Hotplates HD

- \*HD Tops up to 600x450mm\***
- \*Asbestos free insulation\***
- \*Energy control models up to 400°C\***
- \*Thermostat models 50 to 310°C\***
- \*Two year warranty\***



HPR450

FMH Heavy Duty Hotplates conform to the Low Voltage Directive 73/23/EEC and are manufactured to the highest standards with the following features :-

- \* Heavy Duty top \*
- \* Asbestos free insulation \*
- \* Continuously variable heat control by means of an energy regulator \*
- \* Tubular stainless steel sheathed heating element \*
- \* Heater on indicating lamp \*
- \* Robust epoxy coated housing \*
- \* Two year warranty \*
- \* Comprehensive operating manual \*

Special requirements on request.

1) Thermostatic controlled units for more precise temperature control ~ +50 to 310°C and for use with sand bath attachment.

Smaller hotplates with stainless steel tops also available - see Mini Hotplate brochure

MODEL	HPR300	HPR450	HPR600
Temperature range (Surface)	Up to 400°C	Up to 400°C	Up to 400°C
Power supply	230V±10% 50Hz	230V±10% 50Hz	230V±10% 50Hz
Wattage	1500W	2000W	3000W
Temperature control	Energy regulator	Energy regulator	Energy regulator
Nett weight (kg)	12	17	30
Dimensions WxDxH (mm)	370x300x140	520x300x140	670x450x140
Hotplate size (mm)	300x300	450x300	600x450

Specifications :- ambient @ 25°C.

Due to the continuous development of our products, the goods supplied may vary in detail to that illustrated in this brochure.