

Mini hotplate

- *Stainless Tops up to 180x180mm***
- *Spillage lip (MPH180 models only)***
- *Energy control models up to 420°C***
- *Thermostat models 40 to 300°C***
- *Two year warranty***



MHP140

For large heavy duty hotplates
see HPR brochure

FMH Mini Hotplates conform to the Low Voltage Directive 73/23/EEC and are manufactured to the highest standards with the following features :-

- * High grade (304) stainless top *
- * Spillage lip on front (MPH180 models only) *
- * MHP140&180 - Continuously variable heat control by means of an energy regulator *
- * MHP180H - Temperature control by means of a thermostat for more precise temperature control *
- * Tubular stainless steel sheathed heating element *
- * Heater on indicating lamp *
- * Front panel on/off switch *
- * Robust epoxy coated housing *
- * Recessed front control panel *
- * Two year warranty *
- * Comprehensive operating manual *

Options:-

- 1) **Sand bath attachment for MHP180H ONLY**

Special requirements on request.

MODEL	MHP140	MHP180	MHP180H
Temperature range (Surface)	Up to 420°C	Up to 420°C	40 to 300°C
Load Capacity	30 kg	30 kg	30 kg
Power supply	230V±10% 50Hz	230V±10% 50Hz	230V±10% 50Hz
Wattage	500W	840W	840W
Temperature control	Energy regulator	Energy regulator	Hydraulic Thermostat
Nett weight (kg)			
Dimensions WxDxH (mm)	165x210x140	180x224x140	180x240x140
Hotplate size (mm)	140x140	180x180 Spillage lip	180x180 Spillage lip

MHP180H Ideal for use with a Sand bath attachment

Specifications :- ambient @ 25°C.