

# JIPO PORCELAIN

## CRUCIBLES, HIGH FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	19	5	JIP 100301	10
35	22	10	JIP 100351	10
40	25	17	JIP 100401	10
45	28	21	JIP 100451	10
50	32	34	JIP 100501	10
60	38	62	JIP 100601	10
70	44	91	JIP 100701	10



## CRUCIBLES, MIDDLE FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
35	28	12	JIP 200351	10
40	32	20	JIP 200401	10
45	36	30	JIP 200451	10
50	40	45	JIP 200501	10
60	48	80	JIP 200601	10
70	56	120	JIP 200701	10



## CRUCIBLES, LOW FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	38	15	JIP 300301	10
35	44	26	JIP 300351	10
40	50	35	JIP 300401	10
45	56	50	JIP 300451	10
50	62	72	JIP 300501	10
60	75	130	JIP 300601	10



## LIDS FOR CRUCIBLES

Ø (mm) crucible	Ø (mm) lid	Code	Pack of
30	34	JIP 000301	10
35	39	JIP 000351	10
40	44	JIP 000401	10
45	49	JIP 000451	10
50	54	JIP 000501	10
60	64	JIP 000601	10
70	74	JIP 000701	10



## GOOCH CRUCIBLES

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	36	17	JIP 234101	10
35	43	25	JIP 234301	10
39	49	35	JIP 234401	10
60	71	120	JIP 234501	10



## FILTER CRUCIBLES,

with pourous base



H (mm)	Vol (ml)	Porosity	Code	Pack of
28	8	P1 = 6 $\mu$	JIP 400256	10
28	8	P2 = 7 $\mu$	JIP 400257	10
28	8	P2 = 8 $\mu$	JIP 400258	10
35	15	P1 = 6 $\mu$	JIP 400306	10
35	15	P2 = 7 $\mu$	JIP 400307	10
35	15	P2 = 8 $\mu$	JIP 400308	10
40	25	P1 = 6 $\mu$	JIP 400356	10
40	25	P2 = 7 $\mu$	JIP 400357	10
40	25	P2 = 8 $\mu$	JIP 400358	10
43	35	P1 = 6 $\mu$	JIP 400406	10
43	35	P2 = 7 $\mu$	JIP 400407	10
43	35	P2 = 8 $\mu$	JIP 400408	10
50	50	P1 = 6 $\mu$	JIP 400456	10
50	50	P2 = 7 $\mu$	JIP 400457	10
50	50	P2 = 8 $\mu$	JIP 400458	10

## BUECHNER FUNNELS



Nom. size	Stem length	Hole $\varnothing$	Code	Pack of
45	43	1	JIP 237101	10
55	64	1	JIP 237201	10
70	64	1	JIP 237301	8
90	71	2	JIP 237401	6
110	83	2	JIP 237501	1
125	85	2	JIP 237601	1
150	96	2	JIP 237701	1
185	117	2	JIP 237801	1
240	117	2	JIP 237901	1
270	130	2	JIP 237921	1
320	133	2	JIP 237991	1

## DESICCATOR PLATES, ACC. TO CSN.



$\varnothing$ (mm)	for Desic.	Code	Pack of
90	100	JIP 238101	10
140	150	JIP 238331	10
189	200	JIP 238521	10
241	250	JIP 238821	10
290	300	JIP 238991	10

## DESICCATOR PLATES, ACC. TO DIN.



$\varnothing$ (mm)	for Desic.	Code	Pack of
90	100	JIP 138101	10
140	150	JIP 138201	10
190	200	JIP 138301	10
235	250	JIP 138401	10
280	300	JIP 138501	10

## EVAPORATION DISH, ROUND BOTTOM

deep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
42	19	10	JIP 206001	20
51	24	24	JIP 206101	20
61	28	36	JIP 206121	20
77	33	72	JIP 206201	10
81	38	107	JIP 206221	10
94	42	140	JIP 206301	10
103	45	190	JIP 206321	10
112	47	232	JIP 206401	10
132	55	390	JIP 206421	10
138	61	500	JIP 206501	10
148	66	590	JIP 206601	10
170	74	920	JIP 206621	10
180	84	1135	JIP 206701	5
202	86	1500	JIP 206801	5
320	140	6000	JIP 206991	1



## EVAPORATION DISH, ROUND BOTTOM,

semideep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
72	26	41	JIP 207101	20
83	27	62	JIP 207201	20
96	30	110	JIP 207301	20
110	37	160	JIP 207401	10
117	41	215	JIP 207421	10
127	42	240	JIP 207501	10
142	46	325	JIP 207521	10
154	50	420	JIP 207601	10
170	55	558	JIP 207701	5



## EVAPORATION DISH, FLAT BOTTOM,

with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
54	22	22	JIP 274101	20
70	30	62	JIP 274121	20
86	33	93	JIP 274201	20
98	40	154	JIP 274301	10
112	50	265	JIP 274401	10
126	53	330	JIP 274501	10
145	55	450	JIP 274601	10
151	63	620	JIP 274701	10
169	68	800	JIP 274801	10
195	74	1200	JIP 274901	5
254	105	3000	JIP 274991	5



## MORTARS WITH LIP, GLAZED



Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
60	32	25	JIP 211001	10
70	37	30	JIP 211101	10
81	42	50	JIP 211111	5
91	46	75	JIP 211121	5
100	50	110	JIP 211201	5
113	56	150	JIP 211221	5
125	64	220	JIP 211301	5
144	71	325	JIP 211321	1
163	82	500	JIP 211401	1
183	90	730	JIP 211501	1
215	104	1135	JIP 211601	1
257	127	2250	JIP 211701	1
298	150	4300	JIP 211801	1

## PESTLES, GLAZED



Ø (mm)	length (mm)	Code	Pack of
15	52	JIP 213001	10
23	84	JIP 213101	5
26	103	JIP 213201	5
37	132	JIP 213301	5
44	158	JIP 213401	1
58	210	JIP 213501	1
73	255	JIP 213601	1
86	305	JIP 213701	1
103	370	JIP 213801	1

## MORTARS WITH LIP, ROUGH



Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
60	32	25	JIP 211000	10
70	37	30	JIP 211100	10
81	42	50	JIP 211110	5
91	46	75	JIP 211120	5
100	50	110	JIP 211200	5
113	56	150	JIP 211220	5
125	64	220	JIP 211300	5
144	71	325	JIP 211320	1
163	82	500	JIP 211400	1
183	90	730	JIP 211500	1
215	104	1135	JIP 211600	1
257	127	2250	JIP 211700	1
298	150	4300	JIP 211800	1

## PESTLES, ROUGH



Ø (mm)	length (mm)	Code	Pack of
15	52	JIP 213000	10
23	84	JIP 213100	5
26	103	JIP 213200	5
37	132	JIP 213300	5
44	158	JIP 213400	1
58	210	JIP 213500	1
73	255	JIP 213600	1
86	305	JIP 213700	1
103	370	JIP 213800	1